

Decision Chart: CANDOR Product Application



	The Blocking Solution	LowCross-Buffer®	LowCross-Buffer-HRP®	Liquid Plate Sealer®	Antibody Stabilizer	HRP Protector™
Blocking of <i>surface non-specific binding sides</i>	+++	-	-	++	-	-
Blocking of <i>interference in solution/matrix effects</i>	-	+++	+++	-	-	-
Preserving of <i>antibody-enzyme conjugates</i>	-	-	++	-	++	+++
Stabilization of <i>coated proteins/antibodies</i>	+	-	-	+++	-	-
Stabilization of <i>diluted proteins/antibodies</i>	-	+	+++	-	+++	+

KEY: - not recommended + recommended